

## BED & MATTRESS DIMENSIONS

MODEL NAME	MODEL VARIANT	PRODUCT CODE	INTERNAL DIMENSIONS (a)	EXTERNAL DIMENSIONS (b)	MATTRESS DIMENSIONS (c)	SAFE WORKING LOAD (d)
<b>HAMPTON</b>	Standard	CB12TR-200-90	200cm x 90cm	207cm x 101cm	194cm x 83cm x 14.5cm	190kg
	Standard metal rails	CB22SR-200-90	200cm x 90cm	215cm x 93cm	194cm x 83cm x 14.5cm	200kg
<b>HADLEY</b>	Full length metal rails	CB22SRF-200-90	200cm x 90cm	215cm x 93cm	194cm x 83cm x 14.5cm	200kg
	Full length timber rails	CB23TR-200-90	200cm x 90cm	217cm x 101cm	194cm x 83cm x 14.5cm	200kg
<b>PRO-CARE</b>	Standard	CB1-200-90	200cm x 90cm	207cm x 102cm	194cm x 83cm x 14.5cm	190kg
	Wider	CB1-200-120	200cm x 120cm	207cm x 132cm	194cm x 113cm x 14.5cm	190kg
	Longer	CB1-220-90	220cm x 90cm	227cm x 102cm	214cm x 83cm x 14.5cm	190kg
	Shorter	CB1-180-90	180cm x 90cm	187cm x 102cm	174cm x 83cm x 14.5cm	190kg
	Longer & Wider	CB1-220-120	220cm x 120cm	227cm x 132cm	214cm x 113cm x 14.5cm	190kg
<b>PRO-COMM-LOW</b>	Standard	CB2L-200-90	200cm x 90cm	217cm x 101cm	194cm x 83cm x 14.5cm	200kg
	Wider	CB2L-200-120	200cm x 120cm	217cm x 131cm	194cm x 113cm x 14.5cm	200kg
	Longer	CB2L-220-90	220cm x 90cm	237cm x 101cm	214cm x 83cm x 14.5cm	200kg
	Longer & Wider	CB2L-220-120	220cm x 120cm	237cm x 131cm	214cm x 113cm x 14.5cm	200kg
	Double bed	CB2BDuo-200-180	200cm x 180cm	217cm x 198cm	194cm x 83cm x 14.5cm **	400kg
<b>PRO-BARIO</b>	Standard	CB2BL-200-120	200cm x 120cm	217cm x 131cm	194cm x 113cm x 17.8cm	260kg
	Longer	CB2BL-220-120	220cm x 120cm	237cm x 131cm	214cm x 113cm x 17.8cm	260kg
	Max	CB2BLMax-200-120	200cm x 120cm	217cm x 131cm	194cm x 113cm x 17.8cm	325kg
<b>HARVINGTON</b>	Standard	CB3E-200-90	200cm x 90cm	207cm x 109cm	194cm x 88cm x 14.5cm	260kg
	with fixed leg-lift	CB3ST-200-90	200cm x 90cm		194cm x 83cm x 14.5cm	260kg
<b>STANDING BED</b>	with knee-break adjustment	CB3STKB-200-90	200cm x 90cm		194cm x 83cm x 14.5cm	260kg

**Key**

- (a) between siderails and between head & foot boards    (c) actual dimensions    (d) NOT max. user weight    \*\* each mattress  
 (b) assembled bed size. Note: use of Trendelenburg function (where fitted) can increase the external length by up to 24cm